

Pride Jazzy® Elite ES.



Features:

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- Easy-access free-wheel lever
- 60" lap belt

Pride Jazzy® Elite ES specifications:

Weight capacity	300 lbs.
Maximum speed ¹	Up to 4 mph
Length ²	34.75"
Width ²	22.75"
Turning radius	24.75"
Range ^{2,4}	Up to 13.5 miles

Seat sizes:

High-back seat sizes	16" x 16"-18", 18" x 18"-20", 20" x 18"-20"
----------------------	--

Technical:

Drive wheels	9" solid
Rear casters	6" solid casters
Front wheels	3" solid anti-tips
Ground clearance ²	1.5"
Drivetrain	Two-motor, in-line, front-wheel drive
Battery ^{3,5}	(2) 12V, U-1
Battery charger	Off-board, 3.5 amp
Battery weight	25 lbs. each

Colors:

Candy Apple Red	Viper Blue
-----------------	------------

An experience so superior, it can only be elite. The Jazzy® Elite ES features large, front-drive wheels to absorb the bumps and uneven surfaces of everyday life. The comfortable, high-back seat gives you the support you need all day, every day.

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
³ AGM or gel cell type required.
⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
⁵ Battery weight may vary based on manufacturer.

To pre-qualify for a power chair or scooter at little or no cost to you.
 Call toll-free: 866.387.2668, lincare.com