



**Getting to Know Vitera™**

- Introduction 
- Fitting 
- Cleaning 

**Vitera Support Materials**

Sizing Guide



Product Brochure



User Instructions



myMask App



**Mask Parts Replacement Guidelines**

**30 Days**



**Seal**

- Replace every 30 days or when the following occurs:
- Constant leaking
  - Discolored
  - Torn


**3 Months**



**Frame**

- Replace every 3 months or when the following occurs:
- Cracked
  - Discolored

**6 Months**



**Headgear**

- Replace every 6 months or when the following occurs:
- Discolored
  - Hooks or fasteners fall apart
  - Overstretched beyond adjustment