

DeVilbiss Traveler® Portable Compressor Nebulizer System

HCPCS: E0570,A7005,A7013



FIGURE A



FIGURE B



FIGURE C

- Delivers a consistent and concentrated treatment of small particles that penetrate deeply into the lungs, delivering medication where its needed most
- Easy, one-touch operation
- Fast treatment time (approximately 10 minutes) with DeVilbiss 900 Series Reusable Nebulizer encourages patient compliance
- Can be used with any manufacturers nebulizer treatment times will vary
- 45-minute battery (NiMH) life far outlasts treatment time
- Includes DeVilbiss Traveler® compressor, disposable and reusable nebulizers, five extra filters, extra filter door, tubing adapter, adult and pediatric masks, AC adapter, DC car adapter, carrying case, instruction guide and instructional DVD (Figure A)

PRODUCTS	Description
6910D-DR	PORT.COMP. TRAVELER W/O BATT
6910P-DR	Traveler w/ Btty

SPECIFICATIONS	Description
Battery Life	45 Minutes
Battery Warranty	1 Year on Battery
Compressor Warranty	3 Years
Dimensions	3.34" (W) x 2.2" (H) x 3.8" (D) Without Battery 3.34" (W) x 2.2" (H) x 5.3" (D) With Battery
Electrical	100-240V, 50-60Hz, 12VDC, 1.5A
Free Air Flow	6 lpm or Greater
MMAD	5 Microns
Max Pressure	30 psig or Greater
Sound Level	54 dBA
Weight	.8 lbs without battery 1.5 lbs with battery

ACCESSORIES	Description
3655D-621	VixOne Disp. Nebulizer, CS 50, CS
6910D-601	REPLACEMENT BATTERY PACK, EA
DL1050	PEDIATRIC CHARACTER MASK (50-PK), CS
6910D-602	AC ADAPTER, EA
6910D-603	DC ACCESSORY POWER CORD, EA
6910D-605	INLET FILTERS (5) & COVER (1), EA